

Kombinationer Lille Viking 2000/2000+



Rullstol	Modell Lille Viking samt begränsningar	Monteringsats	Länkhjul
	Rullstolar från Etac		
Cross, Cross XL, Twin, Transit SB 37,5–47,5	Lille Viking 2000/2000+ - UB202107/UB202207	UB200624K	175–200 mm
Cross, Cross XL, Twin, Transit SB 50–60	Lille Viking 2000/2000+ - UB202107/UB202207	UB200624L	175–200 mm
Next SB 37,5–47,5	Lille Viking 2000/2000+ - UB202107/UB202207	UB200624N	175–200 mm
Next SB 50–60	Lille Viking 2000/2000+ - UB202107/UB202207	12610	175–200 mm
Next Comfort SB 37,5–50	Lille Viking 2000/2000+ - UB202107/UB202207	UB200624N	175–200 mm
	OBS! På Next Comfort kan drivhjul ej placeras i bakersta balansläget		
M100/M100TR (fr. sitsbredd 40 och sitshöjd 45)	Lille Viking 2000/2000+ - UB202107/UB202207	UB200637A	175–200 mm
	Rullstolar från Invacare		
Rea Focus	Lille Viking 2000/2000+ - UB202107/UB202207	UB200633F	175–200 mm
Rea Spirea	Lille Viking 2000/2000+ - UB202107/UB202207	UB200637A	175–200 mm
Rea Adapt	Lille Viking 2000/2000+ - UB202107/UB202207	UB200637A	175–200 mm
Rea Bellis	Lille Viking 2000/2000+ - UB202107/UB202207	UB200633K	175–200 mm
	Gäller Bellis med 24" drivhjul i högsta sitthöjd		
Rea Assist (gamla modellen)	Lille Viking Flat 2000/2000+ - UB202011/UB445010	UB200637A	175–200 mm
Asalea/Asalea Assist	Lille Viking Flat 2000/2000+ - UB202011/UB445010	UB200629D	175–200 mm
Rea 702	Lille Viking 2000/2000+ - UB202107/UB202207	UB200635A	175–200 mm
Rea 702 TV	Lille Viking 2000/2000+ - UB202107/UB202207	UB200635A	175–200 mm
Rea 703 LX	Lille Viking 2000/2000+ - UB202107/UB202207	UB200635A	175–200 mm
Rea 150	Lille Viking 2000/2000+ - UB202107/UB202207	UB200635A	175–200 mm
	Rullstolar från Hantverksdesign		
HD 500	Lille Viking 2000/2000+ - UB202107/UB202207	UB200626A	175–200 mm
HD 600/650	Lille Viking 2000/2000+ - UB202107/UB202207	UB200626B	175–200 mm
HD 500 Humlan	Lille Viking 2000/2000+ - UB202107/UB202207	UB200626G	175–200 mm
HD 650 Humlan	Lille Viking 2000/2000+ - UB202107/UB202207	UB200626H	175–200 mm

Rullstol	Modell Lille Viking samt begränsningar	Monteringsatts	Länkhjul
	Rullstolar från Handicare		
Impera/Legend/Emineo	Lille Viking 2000/2000+ - UB202107/UB202207	UB200626D	175–200 mm
Cirrus	Lille Viking 2000/2000+ - UB202107/UB202207	UB200626	175–200 mm
Excello	Lille Viking 2000/2000+ - UB202107/UB202207	UB200625K	175–200 mm
Bilanx	Lille Viking 2000/2000+ - UB202107/UB202207	UB200625K	175–200 mm
Exigo	Lille Viking 2000/2000+ - UB202107/UB202207	UB200637A	175–200 mm
	Rullstolar från Alu Rehab		
Netti SB 38–45	Lille Viking 2000/2000+ - UB202107/UB202207	UB200618A	Se monteringsanvisning
Netti SB 48–53	Lille Viking 2000/2000+ - UB202107/UB202207	UB200618	Se monteringsanvisning
Netti 4U SB 40–50	Lille Viking Flat 2000/2000+ - UB202011/UB445010	UB200618A	Se monteringsanvisning
	Rullstolar från Sunrise Medical		
Quickie Classic	Lille Viking 2000/2000+ - UB202107/UB202207	UB200636	175–200 mm
Avantgarde SB 52,5	Lille Viking 2000/2000+ - UB202107/UB202207	UB200637A	175–200 mm
Avantgarde SB från 55	Lille Viking 2000/2000+ - UB202107/UB202207	UB200637B	175–200 mm

74543E

2011-06-27

Etac Sverige AB

Box 203, 334 24 Anderstorp

TELEFON KUNDSERVICE 0371-58 73 00 FAX KUNDSERVICE 0371-58 73 90

E-MAIL info@etac.se HEMSIDA www.etac.se

 etac®